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High-Sensitivity Particle Size Analysis of Colloidal Suspensions (CMP Slurries) by SPOS: Where Less Becomes More

Abstract: The technique of single particle optical sensing (SPOS) offers high resolution and sensitivity in particle size analysis compared with ensemble techniques like Laser Diffraction, gives a number of samples, describes the necessity of using SPOS in CMP process.

Keywords: Single Particle Optical Sensing (SPOS); High resolution; Laser Diffraction; CMP; Particle size analysis

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
Figure 1: a. Population Distributions of Cerium Oxide Slurry 1 (circles) and Slurry 2 (triangles);
b. Volume-Weighted PSDs for Slurry 1 and Slurry 2.

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
Figure 2: Volume-Weighted PSDs of Slurries 1 (circles-unsounded, triangles-sounded) and 2 (diamonds) obtained by Laser Diffraction. a. index of refraction = $1.65 + 0.01i$ b. index of refraction = $1.65 + 0.10i$.

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
Figure 3: Alumina CMP diluted in Sulfuric Acid (circles) and Nitric Acid (squares); a. Number-Weighted PSDs; b. Volume-Weighted PSDs.

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
Figure 4: PSDs of a Silica slurry after being pumped through a recirculation system for 40 hours.

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- [1] D.F. Nicoli, D.C. McKenzie and J.S. Wu, "Application of dynamic light scattering to particle size analysis of macromolecules", *American Laboratory*, Volume 23, No. 17, pp. 32-40 (Nov. 1991)
- [2] D.F. Nicoli, J.S. Wu, Y.J. Chang, D.C. McKenzie and K. Hasapidis, "Automatic, high-resolution particle size analysis by single-particle optical sensing", *American Laboratory*, Volume 24, No. 11, pp. 39-44 (July 1992)
- [3] D.F. Nicoli, K. Hasapidis and D.C. McKenzie, "High Resolution Particle Size Analysis by a Combination of DLS and SPOS", *Lab Asia*, Volume 3, Issue 5, pp. 38-40 (Oct/Nov 1996)
- [4] D.F. Nicoli, J.S. Wu, Y.J. Chang, V. Ovod and K. Hasapidis, "Zeta Potential Analysis of Colloidal Systems by ELS, with Combined DLS Sizing Capability", *International Labmate*, Volume XXII, Issue I, pp. 20-22 (Feb/Mar 1997)
- [5] D.F. Nicoli, J.S. Wu, Y.J. Chang, V. Ovod and K. Hasapidis, "Zeta potential and particle size analysis of colloids using ELS and DLS", *American Laboratory (News Edition)*, Volume 29, No. 19, p. 12 (Sept. 1997)
- [6] D. Nicoli, K. Hasapidis, P. O'Hagan, G. Pokrajac and B. Schade, "Particle size analysis of colloidal suspensions by SPOS compared to DLS: A sensitive indicator of quality and stability", *American Laboratory*, Volume 33, No. 1, pp. 32-39 (Jan. 2001)

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