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The pigments formulated into inkjet inks are typically dispersed to small particle sizes (between about 50 and 200 nm, depending on the application) and need to be made colloiddally stable. The colloiddal stability can be achieved either by surface modification to form an adequate surface charge (zeta potential), or by adsorption of certain compounds on the surface of pigment particles (steric stabilization).

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The size of the pigment particles is critical because large particles may

plug the jets and channels, causing damage to the print head. Controlling the large particle content ($> 0.5\text{-}1.0\ \mu\text{m}$) requires having a technique that is sensitive to a small number of large particles—the tail of the distribution.

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Particle Sizing and Counting Techniques

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There are several methods for determining the mean particle size of ink dispersions such as dynamic light scattering (DLS), but most are not capable of determining small amounts of oversized material. Methods that are based on sizing and counting individual particles are particularly well-suited for this type of analysis, where even small amounts of outliers are sized and counted in the process.

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The AccuSizer single particle optical sizing (SPOS) system is ideal for quantifying the size and concentration of large particle tails present in

inkjet inks. Depending on the sensors incorporated the system can cover a dynamic range of 0.15-400 μm . The system shown in Figure 1 includes the standard LE400 light extinction and scattering sensor that measures from 0.5-400 μm mounted in the AD sampler that provides automated dilution of the sample to the optimum concentration for the measurement. AccuSizer 的单颗粒计数技术(SPOS)可以很好地对油墨中大颗粒定量大小和浓度，具体取决于传感器的检测范围。系统如图 1 所示，包括标准的 LE400 消光和散射传感器——可以测定加入到可自动稀释样液的 AD 采样器的 0.5-400 μm 之间的颗粒，也可以测定试验的最佳浓度。

图 1: AccuSizer 780 AD

Application Example 1: Effect of Stirring

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There are many factors that can impact pigment dispersion one of which

is stirring time. It is necessary to determine an optimal stirring time to reduce the number of oversized particles within a dispersion. It is also important to monitor stirring time since over homogenization can lead to increased particle size.

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Two pigment dispersions, magenta and cyan, were analyzed to monitor the effects of stirring time on oversized particles. The figures 2 and 3 show the results for the magenta sample that was stirred for 50, 70 and 90 minutes. The tail particles decreased from 4×10^6 particles/mL to 2×10^5 particles/mL.

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
Figure 2: Effect of stirring on magenta sample

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
Figure 3: Expanded view of magenta sample

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Concentration of large particles $> 1 \mu\text{m}$:

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70 min: 5×10^5 particles/mL

90 min: 2×10^5 particles/mL

The cyan sample shown in Figure 4 also undergoes a decrease in the number of oversized particles with stirring time. The concentration of particles decreases from 9 million particles/mL to approximately 3 million particles/mL with just 10 additional minutes of stirring.

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
Figure 4: Effect of stirring on cyan sample

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